

READ BEFORE FIRST USE

MAGNET AND FREQUENCY GENERATOR

OPERATING INSTRUCTIONS AND TECHNICAL SPECIFICATIONS

Operating instructions

Magnetic field and frequency device Slow~Wave®

Contents

1. Introduction
2. Scope of delivery
3. Device description
4. Operation and handling
5. Notes regarding water resistance
6. Cleaning
7. Scientific recognition
8. Warranty

1. Introduction

The Slow~Wave® magnetic field and frequency device is a wellness product that generates magnetic field pulses, frequencies, and light. Its applications are similar to natural processes. Please read these operating instructions carefully before using the device. This device cannot replace a healthy lifestyle, meaning adequate sleep, physical outdoor activities, and a healthy diet. The product cannot replace treatment by a doctor or therapist. For a better understanding of the impact on the human body and the effect of pulsating magnetic fields, refer to our product description.

2. Scope of delivery

- Slow~Wave® magnetic field and frequency device
made of wood,
- Induction charger with cable and
charging plug,
- Lanyard and packaging,
- Operating instructions, product description

3. Device description

The Slow~Wave® wellness device is a miniature pulsating magnetic field and frequency device. It has battery-operated electronics that emit magnetic fields, frequencies, and light of various colour spectra.

The enclosure is made of wood (plum, cherry, or birch).

4. Operation and handling

The Slow~Wave® device is shipped charged. The device should be fully charged before it is worn for the first time (charging time approx. 3-4 hours). To charge the device, connect the charging tray using the supplied cable and charging plug (see Fig. 1).

Position the charging tray with the labelled side down. Then place the Slow~Wave® onto the charging tray with its lettering facing the tray. After a short time you will see a blue light (see Fig. 2).

The battery is charged in approx. 3-4 hours, or when the blue light goes off. You can now carry the device (run time approx. 4-6 days).

Once per minute, the display on the device lights up in different colours for approx. 10 s. If the red light flashes rapidly several times in a row (1x, 2x and then 3x), it indicates that charging will soon be required. If the device has not been charged for a long time, there is no light. It will be activated after a short time when the device is recharged.

Position the device with the light outlet approximately in the middle of your breastbone (on your chest). You can adjust the length as you wish. Of course, you can also replace the lanyard with a leather strap or a chain; you are free to decide.

NOTE:

Magnetic field pulses and frequencies are only generated as long as the battery supplies the power electronics with voltage. We therefore recommend charging the device once a week.

5. Note on the water resistance of the device

The housing of the Slow~Wave® is not waterproof. Therefore, the device must not be worn while showering, bathing, or practising water sports. It can be worn when practising other types of sports. Sweat may cause the wood to change its colour slightly, which is normal for natural wood.

6. Cleaning

Clean the device housing from the outside using a soft, clean cloth. Do not use any harsh chemical household cleaners.

7. Scientific recognition

As a manufacturer, we cannot make any statements on healing properties of the Slow~Wave®. Our statements and information about the effects of the device are based on scientific studies

as well as testimonials from users.

We distance ourselves from any statements and promises regarding healing properties, as well as statements on medical efficacy within the meaning of the German Health Services and Products Advertising Act (HWG) and the German Act on Unfair Competition (UWG).

8. Warranty

Provision:

Our warranty period for material defects is 24 months for devices used privately and 12 months for devices used commercially, starting from the delivery date.

Wear parts, consumables, housing surfaces, chains, and eyelets are excluded from warranty claims.

The warranty becomes void:

- in the event of damage resulting from improper operation and improper use of Slow~Wave®, such as severe mechanical stress, impacts, ingress of water, falls;
- if the device has been tampered with or wear parts, consumables, or spare parts that are not original have been used;
- if changes, extensions or repairs to the device are carried out by unauthorised persons;
- if these operating instructions are not observed when the device is in use;
- when securing connecting elements have been loosened, removed, or separated.

Important information

You should drink enough water (about 2 litres / non-carbonated) when using the magnet and frequency generator. If your fluid intake is too low, you may get headaches and any positive effects of the device will be reduced.

At the beginning, there may be short-term reactions such as slight palpitations, increased blood flow to the face, or even slight restlessness. This should not be misunderstood as a side effect. Rather, it is a sign that regulation processes are starting. Those effects should decrease further on. Refer to our booklet, where we describe effects but also exclusion criteria in more detail.

Technical specifications

Battery/accumulator:

- built-in Li-poly battery, 3.7 V, 200 mAh,
- charging via wireless charging, compatible with Qi (WPC v1.2 Qi Industry Standard)
- charging current approx. 40mA -> approx. 5h charging time with an empty battery
- battery life approx. 500 full charge cycles

Display (red LED) when the battery is low:

- red flashing light between colour sequences,
- 1 flash, then 2, etc. up to 5 flashes with increasing discharge,
then the device goes into power-down mode

Magnetic field properties:

- 3 intensities: approx. 800, 1600, 2000 μ T; each frequency change is selected randomly
 - frequencies: 14 frequencies, 8 Hz - 528 Hz, which change at random
- Light runs through different colour spectra once per minute - duration 11 s

Other information:

- weight: approx. 15 g
- dimensions of the wooden body: (height x width x depth): approx. 35 x 35 x 12 mm

The wooden body is a timeless natural product whose dimensions are influenced by basic physical laws such as temperature and humidity. Each tree has its individual features, such as colour, structure and grain.

This natural and rustic character is what makes wood so special.

Each Slow~Wave® is a unique piece.

We hope you enjoy using your new wellness accessory.